



UNIQUE DESIGN FEATURES

CONIC
VASCULAR

Catalogue N°: 22-0040-5 REV.01

MILVUS CONIC
OTW 3.2 F
PTA SEMI COMPLIANT CATHETER



www.conicvascular.com

MILVUS CONIC OTW 3.2 F

PTA SEMI COMPLIANT CATHETER



- Progressive Balloon Diameter
- Deliver central and uniform dilatation
- Short Tubing & Short Marker Bands
- More Flexibility
- Selective use on bifurcated and long lesions
- Conical stent post-dilatation
- Guide Wire 0.014"
- Guide Catheter 5F (0.058")

BALLOON COMPLIANCE TABLE

Pressure (bar)	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm
4	1.84	2.36	2.72	3.34	3.66	4.13	4.80
6	1.92	2.51	2.99	3.49	3.93	4.36	4.91
8	2.00	2.57	3.06	3.65	4.05	4.54	5.02
10	2.04	2.62	3.13	3.72	4.13	4.57	5.10
12	2.09	2.71	3.21	3.81	4.15	4.62	5.18
14	2.18	2.75	3.26	3.94	4.20	4.68	5.28
16	2.23	2.81	3.34	4.06	4.29	4.77	5.36
18	2.28	2.85	3.37	4.13	4.38	4.82	5.42
20	2.31	2.98	3.40	4.21	4.62	-	-
Nominal Ø	8	8	8	8	8	8	8
RBP	18	18	18	18	16	14	14

PRODUCT PART

Product Part Nr.	Balloon Ø	Balloon length
MC-2.5/2.0x30-9	MC-2.5/2.0x30-13	2.5 / 2.0 mm / 30 mm
MC-3.0/2.5x30-9	MC-3.0/2.5x30-13	3.0 / 2.5 mm / 30 mm
MC-3.5/3.0x30-9	MC-3.5/3.0x30-13	3.5 / 3.0 mm / 30 mm
MC-4.0/3.5x30-9	MC-4.0/3.5x30-13	4.0 / 3.5 mm / 30 mm
MC-4.5/4.0x30-9	MC-4.5/4.0x30-13	4.5 / 4.0 mm / 30 mm
MC-5.0/4.5x30-9	MC-5.0/4.5x30-13	5.0 / 4.5 mm / 30 mm

MC-2.5/2.0x40-9	MC-2.5/2.0x40-13	2.5 / 2.0 mm / 40 mm
MC-3.0/2.5x40-9	MC-3.0/2.5x40-13	3.0 / 2.5 mm / 40 mm
MC-3.5/3.0x40-9	MC-3.5/3.0x40-13	3.5 / 3.0 mm / 40 mm
MC-4.0/3.5x40-9	MC-4.0/3.5x40-13	4.0 / 3.5 mm / 40 mm
MC-4.5/4.0x40-9	MC-4.5/4.0x40-13	4.5 / 4.0 mm / 40 mm
MC-5.0/4.5x40-9	MC-5.0/4.5x40-13	5.0 / 4.5 mm / 40 mm

PRODUCT PART

Product Part Nr.	Balloon Ø	Balloon length
MC-4.0/3.0x80-9	MC-4.0/3.0x80-13	4.0 / 3.0 mm / 80 mm
MC-4.5/3.5x80-9	MC-4.5/3.5x80-13	4.5 / 3.5 mm / 80 mm
MC-5.0/4.0x80-9	MC-5.0/4.0x80-13	5.0 / 4.0 mm / 80 mm

MC-4.0/2.5x100-9	MC-4.0/2.5x100-13	4.0 / 2.5 mm / 100mm
MC-4.5/3.0x100-9	MC-4.5/3.0x100-13	4.5 / 3.0 mm / 100mm
MC-5.0/3.5x100-9	MC-5.0/3.5x100-13	5.0 / 3.5 mm / 100mm

MC-4.0/2.0x120-9	MC-4.0/2.0x120-13	4.0 / 2.0 mm / 120mm
MC-4.5/2.5x120-9	MC-4.5/2.5x120-13	4.5 / 2.5 mm / 120mm
MC-5.0/3.0x120-9	MC-5.0/3.0x120-13	5.0 / 3.0 mm / 120mm

